

# BJ TRADE+LINK



## PRODUCT CATALOGUE 2023

# TABLE OF CONTENTS

**Spine Boards & Stretchers** ..... 2



## **SPINE BOARDS & STRETCHERS**

## STRETCHER PVC DOUBLE FOLD INCL. CARRY BAG



PVC Double Fold stretcher with aluminium legs including Carry Bag.

Aluminium alloy folding stretcher used for moving patients who require medical care. Easy storage. Lightweight and strong aluminium alloy frame. Integrated patient restraints to secure patients to the device while transporting.

Product dimensions (L×W×H): 220cm x 55cm x 15cm

Folded dimension: 94cm x 20cm x 9cm

Packing...

[Read More](#)

**SKU:** STRETCHER004

**Price:**





# INDEX

## **S**

Stretcher PVC Double Fold incl. carry bag 3